
Product Datasheet

Anti-FGFR4 / CD334 Reference Antibody (orb1806211)


Description

Anti-FGFR4 / CD334 Reference Antibody is


	hulgG1, and the...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	FGFR4 / CD334
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-FGFR4 / CD334 Reference Antibody (U...



The purity of Anti-FGFR4 / CD334 Referen...



Immobilized human FGFR4 FC at 2 µg/mL C...