



Product Datasheet

Anti-SLC40A1 Reference Antibody (orb1806217)

Description	Anti-SLC40A1 Reference Antibody is expressed from CHO. It has a predicted MW of 145.5 kDa.
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	SLC40A1
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-SLC40A1 Reference Antibody (LY29280...)



The purity of Anti-SLC40A1 Reference Ant...