
Product Datasheet

Anti-TNFRSF17 / BCMA / CD269 Reference Antibody (orb1806259)

Description

Anti-TNFRSF17 / BCMA / CD269 Reference

chain type is hulgG...

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

TNFRSF17 / BCMA / CD269

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

 Maintain refrigerated at 2-8°C for up to 2 weeks.
 For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%


MW

145.5 kDa


Expiration Date

12 months from date of receipt.


kDa NR M R



Anti-TNFRSF17 / BCMA / CD269 Reference A...



The purity of Anti-TNFRSF17 / BCMA / CD2...



U266 cells were stained with Anti-TNFRSF...