

Product Datasheet


Anti-OX2R / CD200R1 Reference Antibody (orb1806261)

Description

Anti-OX2R / CD200R1 Reference Antibody is

kDa NR M R


 Anti-OX2R /
CD200R1
Reference
Antibody
(...

 The purity
of Anti-
OX2R /
CD200R1
Refere...

 Immobilized
human
CD200R His
at 2
µg/mL...

Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	OX2R / CD200R1
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	150 kDa
Expiration Date	12 months from date of receipt.