



Product Datasheet

Anti-TNFRSF9 / 4-1BB / CD137 Reference Antibody (orb1806298)


Description	Anti-TNFRSF9 / 4-1BB / CD137 Reference Antibody is expressed from CHO. The heavy chain type is hulgG...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	TNFRSF9 / 4-1BB / CD137
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145 kDa
Expiration Date	12 months from date of receipt.



Anti-TNFRSF9 /
4-1BB / CD137
Reference A...



The purity of
Anti-TNFRSF9 /
4-1BB / CD1...



Immobilized
humanTNFRSF9
/ 4 1BB /
CD137...