

---


## Product Datasheet

### Anti-IL-6Ra / CD126 Reference Antibody (orb1806335)


**Description**

Anti-IL-6Ra / CD126 Reference Antibody is


kDa NR M R



Anti-IL-6Ra / CD126 Reference Antibody (...



The purity of Anti-IL-6Ra / CD126 Refere...



Immobilized human IL-6R His at 2 µg/mL ...

146.32 kDa.

**Reactivity**

Human

**Endotoxins**

&lt; 1 EU/mg

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, FA, FACS, Kinetics

**Target**

IL-6Ra / CD126

**Preservatives**

100 mM Pro-Ac 20 mM Arg pH 5.0

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Monoclonal

**Purity**

&gt; 95%

**MW**

146.32 kDa

**Expiration Date**

12 months from date of receipt.