

---

## Product Datasheet


### Anti-FGFR2 / CD332 Reference Antibody (orb1806349)

**Description**


Anti-FGFR2 / CD332 Reference Antibody is

	hulgG1, and the...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/mg
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FA, FACS, Kinetics
<b>Target</b>	FGFR2 / CD332
<b>Preservatives</b>	100 mM Pro-Ac 20 mM Arg pH 5.0
<b>Storage</b>	-80°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95%
<b>MW</b>	78.2 kDa
<b>Expiration Date</b>	12 months from date of receipt.


kDa NR M R



Anti-FGFR2 / CD332 Reference Antibody (a...



The purity of Anti-FGFR2 / CD332 Referen...



Human FGFR2-β-IIIc (A13) CHOS cells wer...