

---

## Product Datasheet

### Anti-FGFR2 / CD332 Reference Antibody (orb1806350)

**Description**

Anti-FGFR2 / CD332 Reference Antibody is


kDa NR M R



Anti-FGFR2 / CD332 Reference Antibody (b...



The purity of Anti-FGFR2 / CD332 Referen...



SNU-16 cells were stained with Anti-FGFR...

hulgG1, and the...

**Reactivity**

Human

**Endotoxins**

&lt; 1 EU/mg

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, FA, FACS, Kinetics

**Target**

FGFR2 / CD332

**Preservatives**

100 mM Pro-Ac 20 mM Arg pH 5.0

**Storage**

-80°C

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Monoclonal

**Purity**

&gt; 95%

**MW**

145.24 kDa

**Expiration Date**

12 months from date of receipt.