
Product Datasheet


Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody (orb1806356)

Description

Anti-TNFRSF10B / TRAILR2 / CD262 Reference

kDa NR M R

chain type is h...


Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

TNFRSF10B / TRAILR2 / CD262

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

-80°C

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%

MW

144.96 kDa


Expiration Date

12 months from date of receipt.

Anti-TNFRSF10B / TRAILR2 / CD262 Referen...



The purity of Anti-TNFRSF10B / TRAILR2 / ...



Immobilized Cyno DR5 His at 2 µg/mL can...