


Product Datasheet


Anti-DLL3 Reference Antibody (orb1806357)

Description	Anti-DLL3 Reference Antibody is expressed from CHO. It has a predicted MW of 144.16 kDa.
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	DLL3
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	144.16 kDa
Expiration Date	12 months from date of receipt.



Anti-DLL3 Reference Antibody (Dragonfly

...



The purity of Anti-DLL3 Reference Antibo...