

Product Datasheet

Anti-CX3CL1 / Fractalkine Reference Antibody (orb1806362)

Description Anti-CX3CL1 / Fractalkine Reference Antibody is expressed from CHO. The heavy chain type is hulgG2, ...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target CX3CL1 / Fractalkine

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only


Isotype IgG

Clonality Monoclonal

Purity > 95%

MW 150 kDa

Expiration Date 12 months from date of receipt.



Anti-CX3CL1 / Fractalkine Reference Antibody



The purity of Anti-CX3CL1 / Fractalkine

...