
Product Datasheet


Anti-CCN2 / CTGF Reference Antibody (orb1806364)

Description


Anti-CCN2 / CTGF Reference Antibody is expressed

	I...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	CCN2 / CTGF
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	150 kDa
Expiration Date	12 months from date of receipt.


kDa NR M R



Anti-CCN2 / CTGF Reference Antibody (pam...



The purity of Anti-CCN2 / CTGF Reference...



Immobilized human CTGF His at 2 µg/mL c...