

Product Datasheet

Anti-H3L Antibody (orb1806387)

Description	Anti-H3L Antibody synthesised from CHO. The antibody type is VHH-Fc(hulgG1). It has a predicted MW of 146.12 ...
Reactivity	Virus
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	H3L
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	VHH
Clonality	Monoclonal
Purity	> 95%
MW	146.12 kDa
Expiration Date	12 months from date of receipt.