

Product Datasheet

Anti-E8L Antibody (orb1806388)

Description	Anti-E8L Antibody synthesised from CHO. The heavy chain type is hulgG1, and the light chain type is hukappa....
Reactivity	Virus
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	E8L
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	78.64 kDa
Expiration Date	12 months from date of receipt.