

Product Datasheet

Anti-Spike S1 (Omicron BA.2 / BA.2.75 / BA.4 / BA.5) Antibody (orb1806390)

Description	Anti-Spike S1 (Omicron BA.2 / BA.2.75 / BA.4 / BA.5) Antibody is expressed from CHO. The heavy chain...
Reactivity	Virus
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	Spike S1
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.6 kDa
Expiration Date	12 months from date of receipt.