


Product Datasheet

CD315 Antibody (PE) (orb1806409)

Description	CD315, also known as prostaglandin F2 receptor negative regulator, is an approximately 135 kDa trans...
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	HeLa cells
Target	CD315
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	Q9P2B2
Entrez	5738
Expiration Date	12 months from date of receipt.



Separation of A431 cells stained using a...