



Product Datasheet

IL-4 Antibody (orb1806424)


Description	IL-4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Residues 25-153
Target	IL-4
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P05112
Dilution Range	WB: 1:1000, IP: 20 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human IL-4 by western blot....



Detection of human IL-4 by western blot....



Detection of human IL-4 by western blot ...