



Product Datasheet

PVRL1/NECTIN1 Antibody (orb1806450)

Description	PVRL1/NECTIN1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 225 and 275
Target	PVRL1/NECTIN1
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	500 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	Q15223
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Detection of human PVRL1/NECTIN1 by west...



Detection of mouse PVRL1/NECTIN1 by west...