
Product Datasheet

**Biotinylated Human Galectin-9 / LGALS9 Protein,
His, (SPR verified) (orb1806470)**

Description

Biotinylated Human Galectin-9,

293 cells (HEK293). It contains AA Ala 2...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

Galectin-9

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

N-10×His & N-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus, followed by an Avi tag (). The protein has a calculated MW of 39.3 kDa. The protein migrates as 52-58 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Biotinylated Human Galectin-9, His, is expressed from human 293 cells (HEK293). It contains AA Ala 2 - Thr 323 (Accession # O00182-2).

MW

39.3 kDa

Uniprot ID**O00182****Expiration Date**

6 months from date of receipt.