
Product Datasheet

**Rhesus macaque TLR3 / CD283 Protein, Fc Tag
(MALS verified) (orb1806474)**

Description

Rhesus macaque TLR3 Protein, Fc

	293 cells (HEK293). It contains AA Ser 2...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TLR3
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 103.4 kDa. The protein migrates as 130-160 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Rhesus macaque TLR3 Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Ser 23 - Glu 703 (Accession # F7DWP4-1).
MW	103.4 kDa
Uniprot ID	F7DWP4
Expiration Date	6 months from date of receipt.