

Product Datasheet

IgG control Antibody (APC) (orb1806498)

Description	Mouse IgG2b Antibody Isotype Control is a monoclonal antibody recombinantly expressed from human 293...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	APC
Target	IgG control
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 12 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower. Please protect from light and avoid repeated freeze-thaw cycles
Note	For research use only
Application notes	Flow Cytometry.
Isotype	Mouse IgG2b Mouse kappa
Purity	90%
Expiration Date	12 months from date of receipt.