


Product Datasheet

Human CRLF2/TSLP R ELISA Kit (orb180654)

Description	Human CRLF2/TSLP R ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G25-K231
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 62.5 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	42013 MW
Uniprot ID	Q9HC73
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human CRLF2/TSLP R PicoKine ELISA Kit st...