


Product Datasheet

Rat IGFBP5 ELISA Kit (orb180662)

Description	Rat IGFBP5 ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L20-E271
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:100, concentration in serum and plasma is around 180 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	30298 MW
Uniprot ID	P24594
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Rat IGFBP5
PicoKine
ELISA Kit
standard c...