

## Product Datasheet

# Human CD169/SIGLEC-1/Sialoadhesin ELISA Kit (orb180668)

**Description** Human CD169/SIGLEC-1/Sialoadhesin ELISA Kit

**Reactivity** Human

**Range** 156 pg/ml - 10,000 pg/ml

**Tested Applications** ELISA

**Immunogen** Expression system for standard: NS0; Immunogen sequence: S20-Q1641

**Storage** Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

**Note** For research use only

**Application notes** Dilution ratio of 1:1, concentration in serum and plasma is around 4 ng/ml


**Sample Types** cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).

**MW** 182624 MW

**Uniprot ID** [Q9BZZ2](#)

**Sensitivity** 10 pg/ml

**Expiration Date** Please enquire.



Human CD169/SIGLEC-1 PicoKine ELISA Kit ...