


Product Datasheet

Human TGFBR2 / Tgf Beta Receptor II / TGF beta R2 ELISA Kit (orb180671)

Description	Human TGFBR2 / Tgf Beta Receptor II / TGF beta R2 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: I24-D159
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:4, concentration in serum and plasma is around 600 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	64568 MW
Uniprot ID	P37173
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human TGFBR2 PicoKine ELISA Kit standard...