


Product Datasheet

Human CRP/C Reactive Protein ELISA Kit (orb180677)

Description	Human CRP/C Reactive Protein ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F17-P224
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:5000, concentration in serum and plasma is around 10 ug/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	25039 MW
Uniprot ID	P02741
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human CRP PicoKine ELISA Kit standard cu...