

Product Datasheet

Rat CRH(Corticotropin Releasing Hormone) ELISA Kit (orb1806844)

Description Rat CRH(Corticotropin Releasing Hormone) ELISA Kit

Reactivity Rat

Range 0.31-20ng/mL

Target CRH

Note For research use only

Application notes This ELISA kit uses the Competitive-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with Rat CRH. During the reaction, Rat CRH in the sample or standard competes with a fixed amount of Rat CRH on the solid phase supporter for sites on the Biotinylated Detection Ab specific to Rat CRH. Excess conjugate and unbound sample or standard are washed away, and Avidin-Horseradish Peroxidase (HRP) conjugate are added to each micro plate well and incubated. Then a TMB substrate solution is added to each well. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns from blue to yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm $\pm$ 2 nm. The concentration of Rat CRH in tested samples can be calculated by comparing the OD of the samples to the standard curve.

Sample Types serum, plasma, Tissue homogenate and Other biological samples

Sensitivity 0.19 ng/mL

Expiration Date Please enquire.

