

Product Datasheet

Recombinant Bovine Pregnancy-associated glycoprotein 1 (HEK293) (orb1808040)

Description	Recombinant Bovine Pregnancy-associated glycoprotein 1 (HEK293)
Endotoxins	Not analyzed
Target	PAG 1
Form/Appearance	Lyophilized or liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
Uniprot ID	Q29432
Expiration Date	6 months from date of receipt.