

Product Datasheet

Recombinant mouse PD-L1 protein, C-hFc (HEK293) (orb1808046)

Description	Recombinant mouse PD-L1 protein, C-hFc (HEK293)
Endotoxins	Not analyzed
Target	CD274
Form/Appearance	Lyophilized or liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-hFc
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	19-238/290
Purity	> 90% as determined by SDS-PAGE
MW	51.5 kDa
Uniprot ID	Q9NZQ7
Expiration Date	6 months from date of receipt.