

Product Datasheet

Recombinant mouse CD3E protein, C-mFc (HEK293) (orb1808061)

Description	Recombinant mouse CD3E protein, C-mFc (HEK293)
Target	CD3E
Form/Appearance	Lyophilized or liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-mFc
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	23-108/189
MW	37.2 kDa
Uniprot ID	P07766
Expiration Date	6 months from date of receipt.