

## Product Datasheet

# Recombinant human Insulin protein, N-Trx-His (orb1808065)

<b>Description</b>	Recombinant human Insulin protein, N-Trx-His
<b>Endotoxins</b>	Not analyzed
<b>Target</b>	INS
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	> 1mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	N-Trx-His
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Protein Sequence</b>	25-110/110
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>MW</b>	27.2 kDa
<b>CAS Number</b>	11061-68-0
<b>Uniprot ID</b>	<b>P01308</b>
<b>Expiration Date</b>	6 months from date of receipt.