

Product Datasheet

Recombinant human CD56 protein, C-His (HEK293) (orb1808069)

Description	Recombinant human CD56 protein, C-His (HEK293)
Target	NCAM1
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	21-603/761
Purity	> 90% as determined by SDS-PAGE
MW	65.7 kDa
Uniprot ID	P13591
Expiration Date	6 months from date of receipt.