

## Product Datasheet

# Recombinant human EphA3 protein, C-His (HEK293) (orb1808088)

<b>Description</b>	Recombinant human EphA3 protein, C-His (HEK293)
<b>Target</b>	EPHA3
<b>Form/Appearance</b>	Lyophilized or liquid
<b>Concentration</b>	> 0.5mg/ml
<b>Storage</b>	> 90% as determined by SDS-PAGE
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Protein Sequence</b>	21-541/983
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>MW</b>	60 kDa
<b>Uniprot ID</b>	<a href="#">P29320</a>
<b>Expiration Date</b>	6 months from date of receipt.