

Product Datasheet

Recombinant human SNAI2 protein, N-His (orb1808111)

Description	Recombinant human SNAI2 protein, N-His
Endotoxins	Not tested
Target	SNAI2
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	N-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	1-268/268
Purity	> 90% as determined by SDS-PAGE
MW	32.1 kDa
Uniprot ID	O43623
Expiration Date	6 months from date of receipt.