

Product Datasheet

Recombinant mouse Angiopoietin 2 protein, N-His (orb1808135)

Description	Recombinant mouse Angiopoietin 2 protein, N-His
Endotoxins	Not analyzed
Target	ANGPT2
Form/Appearance	Lyophilized or liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	N-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	263-496/496
Purity	> 90% as determined by SDS-PAGE
MW	28.7 kDa
Uniprot ID	O15123
Expiration Date	6 months from date of receipt.