

Product Datasheet

Recombinant human Osteopontin protein, C-His (HEK293) (orb1808138)

Description	Recombinant human Osteopontin protein, C-His (HEK293)
Target	SPP1
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Purity	> 90% as determined by SDS-PAGE
Uniprot ID	P10923
Expiration Date	6 months from date of receipt.