

Product Datasheet

Recombinant mouse GM-CSFR alpha protein, His (HEK293) (orb1808139)

Description	Recombinant mouse GM-CSFR alpha protein, His (HEK293)
Target	CSF2RA
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	30-327/388
Purity	> 90% as determined by SDS-PAGE
MW	32.8 kDa
Uniprot ID	P15509
Expiration Date	6 months from date of receipt.