

Product Datasheet

Recombinant human Podoplanin protein, C-hFc (HEK293) (orb1808148)

Description	Recombinant human Podoplanin protein, C-hFc (HEK293)
Endotoxins	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Target	PDPN
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-hFc
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
MW	38 kDa
Uniprot ID	Q86YL7
Expiration Date	6 months from date of receipt.