

Product Datasheet

bs-0295P-RBITC Antibody (orb1808170)

| | |
|------------------------|---|
| Description | bs-0295P-RBITC Antibody |
| Target | bs-0295P-RBITC |
| Preservatives | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% glycerol. |
| Storage | Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |