

## Product Datasheet

### Mouse ANXA5 ELISA Kit (orb1808371)

|                        |                                           |
|------------------------|-------------------------------------------|
| <b>Description</b>     | Mouse ANXA5 ELISA Kit                     |
| <b>Reactivity</b>      | Mouse                                     |
| <b>Range</b>           | 15.6-1000pg/mL                            |
| <b>Storage</b>         | Store at 4°C & -20°C for 12 months.       |
| <b>Note</b>            | For research use only                     |
| <b>Sample Types</b>    | Serum, plasma or other biological fluids. |
| <b>Sensitivity</b>     | 5.8pg/mL                                  |
| <b>Expiration Date</b> | Please enquire.                           |