

## Product Datasheet

### SB242235 (orb180841)

<b>Description</b>	SB-242235 is a potent and selective p38 MAP kinase inhibitor with IC50 of 1.0 uM.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>19</sub> H <sub>20</sub> F <sub>5</sub> O
<b>MW</b>	353.39
<b>CAS Number</b>	[193746-75-7]
<b>SMILES</b>	<chem>COC1=NC=CC(=N1)C2=C(N=CN2C3CCNCC3)C4=CC=C(C=C4)F</chem>
<b>Expiration Date</b>	12 months from date of receipt.