

## Product Datasheet

### Human PIK3C2a ELISA Kit (orb1808509)

<b>Description</b>	Human PIK3C2a ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3.12-200U/L
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Tissue homogenates or other biological fluids.
<b>Sensitivity</b>	1.6U/L
<b>Expiration Date</b>	Please enquire.