

## Product Datasheet

### Human KLK5 ELISA Kit (orb1808624)

<b>Description</b>	Human KLK5 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78.1-5000pg/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, saliva, tissue homogenates or other biological fluids.
<b>Sensitivity</b>	29pg/mL
<b>Expiration Date</b>	Please enquire.