

## Product Datasheet

### Human HK1 ELISA Kit (orb1808657)

<b>Description</b>	Human HK1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5-4000pg/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates or other biological fluids.
<b>Sensitivity</b>	27.2pg/mL
<b>Expiration Date</b>	Please enquire.