

Product Datasheet

Human GHRL ELISA Kit (orb1808683)

| | |
|------------------------|--|
| Description | Human GHRL ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Storage | Store at 4°C & -20°C for 12 months. |
| Note | For research use only |
| Sample Types | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sensitivity | 0.047ng/mL |
| Expiration Date | Please enquire. |