

## Product Datasheet

### VS-5584 (orb180870)

<b>Description</b>	VS-5584 (SB2343) is a potent and selective dual PI3K/mTOR inhibitor for mTOR, PI3K $\alpha$ / $\beta$ / $\delta$ / $\gamma$ with I...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>17</sub> H <sub>22</sub> N <sub>8</sub> O
<b>MW</b>	354.41
<b>CAS Number</b>	[1246560-33-7]
<b>SMILES</b>	CC1=NC2=C(N=C(N=C2N1C(C)C)N3CCOCC3)C4=CN=C(N=C4)N
<b>Expiration Date</b>	12 months from date of receipt.