

## Product Datasheet

### Human DPP4 ELISA Kit (orb1808730)

<b>Description</b>	Human DPP4 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.2-2000pg/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates or other biological fluids.
<b>Sensitivity</b>	14.8pg/mL
<b>Expiration Date</b>	Please enquire.