

Product Datasheet

Human DKK1 ELISA Kit (orb1808732)

Description	Human DKK1 ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Storage	Store at 4°C & -20°C for 12 months.
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	30pg/mL
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)