

## Product Datasheet

### Human CASP14 ELISA Kit (orb1808797)

|                        |                                                              |
|------------------------|--------------------------------------------------------------|
| <b>Description</b>     | Human CASP14 ELISA Kit                                       |
| <b>Reactivity</b>      | Human                                                        |
| <b>Range</b>           | 0.312-20ng/mL                                                |
| <b>Storage</b>         | Store at 4°C & -20°C for 12 months.                          |
| <b>Note</b>            | For research use only                                        |
| <b>Sample Types</b>    | Tissue homogenates, cell lysates or other biological fluids. |
| <b>Sensitivity</b>     | 0.106ng/mL                                                   |
| <b>Expiration Date</b> | Please enquire.                                              |