

## Product Datasheet

### General DHVD3 ELISA Kit (orb1808932)

<b>Description</b>	General DHVD3 ELISA Kit
<b>Reactivity</b>	All
<b>Range</b>	15.6-1000pg/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma or other biological fluids.
<b>Sensitivity</b>	4.1pg/mL
<b>Expiration Date</b>	Please enquire.